

Clonal Technologies

For more information
contact Clonal Technologies
on 1300 795 335


See our range online at www.clonaltech.com.au

Email: john.dubois@clonaltech.com.au

33 Fulcrum Street
Richlands Queensland 4077

Instant-View® Multi-Panel Urine Test Drug Screen Cassettes

Instant-View® Multi-Panel 6
Test Drug Screen Cassettes
detect DoAs in urine to the
following standards:
Amphetamine (300 ng/ml),
Benzodiazepine (200 ng/ml),
Cocaine (300 ng/ml),
Methamphetamine (300 ng/ml),
MDMA/XTC (300 ng/ml),
Morphine (300 ng/ml), and
THC (50 ng/ml).



Features

- Independently validated to meet Australian Standards 4308
- Quick results (45 secs)
- Very strong bands
- Individual test strips with Test and Control line for each drug class
- Colour-coded strips
- Dip or pipette method
- Sealed cassette and sample well
- Proven history in Australia
- Extensive technical support

Benefits

- Consistent with standards for cut off levels
- High throughput
- Easy interpretation of results
- Easier interpretation and recognition of drug class
- Easier recognition of drug class
- Ease of use, flexibility
- Reduces urine exposure, cleaner, and allows photocopying
- Reliability and confidence
- Full company expertise at your disposal