



Clonal Technologies

For more information
contact Clonal Technologies
on 1300 795 335

See our range online at www.clonaltech.com.au

Email: john.dubois@clonaltech.com.au

33 Fulcrum Street
Richlands Queensland 4077

Vision Test Cup

Vision Test Cups detect DoAs in urine to the following standards:
Amphetamine (300 ng/ml),
Benzodiazepine (200 ng/ml),
Cocaine (300 ng/ml),
Methamphetamine (300 ng/ml),
MDMA/XTC (300 ng/ml),
Morphine (300 ng/ml),
and THC (50 ng/ml).



Features

Independently validated to meet
Australian Standards 4308.

Quick results.

Very strong bands.

Individual test strips with Test and
Control line for each drug class.

Self-sealing urine cup.

In-built sampling tube.

Covered sampling hole in lid.

Temperature strip.

Proven history in Australia.

Extensive technical support.

Benefits

Consistent with standards for cut off levels.

High throughput.

Easy interpretation of results.

Easier interpretation and recognition
of drug class.

Cup is test well; less steps for sampling;
Cleaner; No-mess transport to labs.

No tipping or button push required.

Convenient resampling for laboratory

Prevention of adulteration.

Reliability and confidence.

Full company expertise at your disposal.