



For more information
contact Clonal Technologies
on 1300 795 335

See our range online at www.clonaltech.com.au

Email: john.dubois@clonaltech.com.au

33 Fulcrum Street
Richlands Queensland 4077

Integrated E-Z Split Cup

Urine Test Drug Screen

Integrated E-Z Split Cup Urine Test Drug Screens detect DoAs in urine to the following standards:

Amphetamine (300 ng/ml),
Benzodiazepine (200 ng/ml),
Cocaine (300 ng/ml),
Methamphetamine (300 ng/ml),
Morphine (300 ng/ml),
and THC (50 ng/ml).



Features

Six separate adulteration strips which includes oxidants/PCC, specific gravity, pH, nitrite, glutaraldehyde and creatinine.

Hinged cup lid

Independently validated to meet Australian Standards 4308.

Quick results.

Very strong bands.

Individual test strips with Test and Control line for each drug class.

Urine cup is test well.

Cup lid with split key.

Key activation.

Temperature strip.

Proven history in Australia.

Extensive technical support.

Benefits

Prevention of adulteration.

Lid stays attached.

Consistent with standards for cut off levels.

High throughput.

Easy interpretation of results.

Easier interpretation and recognition of drug class.

Less steps for sampling.

Reduces urine exposure, cleaner. Sample can be saved for labs.

Prevention of adulteration.

Reliability and confidence.

Full company expertise at your disposal.