



Clonal Technologies


Clonal Technologies Pty Ltd ABN 41 104 918 434

33 Fulcrum Street,
Richlands, Qld 4077

Tel 07 3859 2950
Fax 07 3859 2998


Web www.clonaltech.com.au
Email clonal@clonaltech.com.au

Integrated E-Z Split Cup II Urine Test Drug Screen



Integrated E-Z Split Cup II Urine Test Drug Screens detect DoAs in urine to the following standards:

- Amphetamine (300 ng/ml),
- Benzodiazepine (200 ng/ml),
- Cocaine (300 ng/ml),
- Methamphetamine (300 ng/ml),
- Morphine (300 ng/ml),
- and THC (50 ng/ml).



Features

Six separate adulteration strips which includes oxidants/PCC, specific gravity, pH, nitrite, glutaraldehyde and creatinine.

Hinged cup lid

Independently validated to meet Australian Standards 4308.

Quick results.

Very strong bands.

Large flat Test Panel.

Urine cup is test well.

Cup lid with split key.

Key activation.

Temperature strip.

Proven history in Australia.

Extensive technical support.

Benefits

Prevention of adulteration.

Lid stays attached.

Consistent with standards for cut off levels.

High throughput.

Easy interpretation of results.

Results are easily scanned or photocopied.

Less steps for sampling.

Reduces urine exposure, cleaner. Sample can be saved for labs.

Prevention of adulteration.

Reliability and confidence.

Full company expertise at your disposal.